



Patient Care Products

A World of Medical Needs

INTRODUCING FLEXYFOOT

Flexyfoot is a premium dynamic rubber tip designed for confidence and comfort when using a walking aid. Patented air sprung suspension provides cushioning and allows full contact of the rubber tip with the ground, improving grip and stability on all terrains, including slippery and uneven surfaces.



Great on the ground:

Flexyfoot will always maintain full contact with the floor, no matter what angle the walking aid is used at, providing improved grip and stability on all terrains, and in particular with slippery and uneven surfaces. The Flexyfoot rubber tread can also independently rotate through 360°, so full ground contact is still maintained even when twisting or turning around.

Gentle on the joints:

Flexyfoot's patented design works like the suspension system on a four wheel drive. The flexible bellows absorb shock each time the walking aid contacts the ground so that forces transmitted to the upper body through the aid are always cushioned.

Easy to fit:

Flexyfoot looks a lot smarter than a traditional ferrule and is much easier to monitor for wear and replace when necessary. When the rubber tread is worn down to the wear markers (similar to those of a car tyre), the bellows and integral rubber foot are easily screwed off from the Flexyfoot collar and replaced.

Flexyfoot comes in four collar sizes to fit the most common walking aids. It is also perfect for improving the safety of shower chairs, which are used in wet and cambered shower recesses.

19 mm is suitable for most walking sticks, including adjustable and folding walking sticks.

22 mm is suitable for most adjustable crutches.

25 mm is suitable for most walking frames and shower chairs.

16 mm is suitable for most canes and small diameter walking sticks.



The problem...

The Flexyfoot solution...

Why not ask about getting Flexyfoot fitted here today?

australia flexyfoot

Walk with confidence!